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IFU
Rongeurs

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
Contents

IFU – Instructions for use Open Up®-System RW 5xxx O	30
IFU – Instructions for use Open Up®-System RW 49xx O	33
IFU – Instructions for use Open Up®-System RW 3xxx O / 46xx O / RW 47xx O / RW 48xx O	36
1 Scope	39
2 Intended use and indication	41
2.1 Intended use	41
2.2 Indication(s)	41
2.3 Contraindication(s)	42
3 Safety	42
3.1 Information about warnings	42
3.1.1 Warning symbols	42
3.2 General safety instructions	42
3.3 Checking before use	43
3.4 Safety during use	43
3.5 Reportable incidents	46
4 Reprocessing procedure	46
4.1 General instructions	46
4.2 Maximum number of reprocessing cycles	46
4.3 Initial treatment at the site of use	46
4.4 Precleaning	47
4.5 Automated cleaning and disinfection	48
4.6 Inspection	49
4.6.1 Visual inspection	49
4.6.2 Functional inspection	50
4.7 Packaging	50
4.8 Steam sterilization	51
5 Storage	51
5.1 Sterile, packaged instruments	51
5.2 Nonsterile instruments	51
6 Technical service	52
6.1 Manufacturer/ Service address:	52
7 Disposal	53
8 Symbols	54
9 Trademarks	55
9.1 Trademarks of Raimund Wenzler GmbH	55
9.2 Trademarks of third-party manufacturers	55

IFU – Instructions for use

Open Up[®]-System RW 5xxx O


Scope

Article groups

REF

RW 5xxx O


Open Up[®]-System

- 1** Footplate
- 2** Main shaft
- 3** Sliding shaft
- 4** Push button
- 5** Handle part, movable
- 6** Handle part, fixed
- 7** External guide with ejector
- 8** Arrows


Opening up for the cleaning position

- ▶ Press and hold down the push button (4)
- ▶ Pull the movable handle part (5) down to the stop position

- ▶ Pull the sliding shaft (3) slightly back and then up

- ▶ The instrument is in the opened-out cleaning position


Locking

- ▶ Insert the sliding shaft (3) into the external guide (7) of the main shaft

- ▶ Push the sliding shaft (3) forward until the arrows are on top of one another

- ▶ Move the movable handle part (5) up until the locking system audibly engages.
- ▶ The instrument is ready for use.

IFU – Instructions for use Open Up®-System RW 49xx O


Scope:

Article groups

REF

RW 49xx O


Open Up®-System

- 1** Jaw section (movable and fixed)
- 2** Main shaft
- 3** Sliding shaft
- 4** Push button
- 5** Handle part, movable
- 6** Handle part, fixed
- 7** Guide (inside sliding shaft)
- 8** Arrows

Opening up for the cleaning position


- ▶ Press and hold down the push button (4)
- ▶ Pull the movable handle part (5) down to the stop position

- ▶ Pull the sliding shaft (3) slightly back and then up
- ▶ The instrument is in the opened-out cleaning position

Locking

- ▶ Insert the sliding shaft (3) into the internal guide (7) of the main shaft

- ▶ Push the sliding shaft (3) forward until the arrows are on top of one another

- ▶ Move the movable handle part (5) up until the locking system audibly engages.
- ▶ The instrument is ready for use.

IFU – Instructions for use


Open Up®-System RW 3xxx O / 46xx O / RW 47xx O / RW 48xx O

Scope:

Article groups

REF

RW 3xxx O / RW 46xx O / RW 47xx O / RW 48xx O
Open Up®-System


- 1 Jaw section (movable and fixed)
- 2 Main shaft
- 3 Sliding shaft
- 4 Push button
- 5 Handle part, movable
- 6 Handle part, fixed
- 7 Guide (inside sliding shaft)
- 8 Arrows

Opening up for the cleaning position

- ▶ Bring the movable handle part (5) into the opened position (jaw section is open)


- ▶ Press and hold down the push button (4)
- ▶ Pull the movable handle part (5) down to the stop position
The push button (4) can be released


- ▶ Pull the sliding shaft (3) slightly back and then up
- ▶ The instrument is in the opened-out cleaning position

Locking

- ▶ Insert the sliding shaft (3) into the internal guide of the main shaft (7)

- ▶ Place both arrows (8) on top of one another and then move the handle part (5) up until the locking system audibly engages.
- ▶ The instrument is ready for use.

1 Scope

These instructions for use apply to the following article groups

Catalogue number REF		Hard bone	Porous bone structures/ cartilage	Soft tissue
RW 208X RW 209X	Decker rongeur		x	x
RW 36XX RW 37XX	Caspar rongeur		x	x
RW 471X RW 481X RW 452X RW 462X RW 472X RW 482X	Cushing IVD rongeur ...		x	x
RW 453X RW 463X RW 473X RW 483X	Love-Grünwald IVD rongeur ...		x	x
RW 454X RW 464X RW 474X RW 484X RW 455X RW 465X RW 475X	Spurling IVD rongeur ...		x	x
RW 480X	Schlesinger rongeur		x	x
RW 487X RW 488X	Peapod rongeur		x	x
RW 4848	Hoen rongeur		x	x
RW 4849	Williams rongeur		x	x
RW 4970	Oldberg rongeur		x	x
RW 49XX	Ferris-Smith rongeur ...		x	x
RW 50XX RW 51XX RW 52XX RW 53XX RW 57XX RW 58XX RW 59XX	Ferris-Smith Kerrison rongeur ...	Limitations on use: Only with footplate version: regular: 	No limitations on use: Footplate version: regular and thin (DFP):  	No limitations on use: Footplate version: regular and thin (DFP):  

Catalogue number REF			Hard bone	Porous bone structures / cartilage	Soft tissue
RW 52XX	RW 56XX	Ferris–Smith Kerrison bayonet rongeur	Limitations on use: Only with footplate version: regular: 	No limitations on use: Footplate version: regular and thin (DFP):  	No limitations on use: Footplate version: regular and thin (DFP):  
RW 893X	RW 894X	Leksell rongeur	x	x	x
RW 8950		Beyer rongeur	x	x	x

The rongeurs are available with different variable dimensions:

- ▶ Shaft length
- ▶ Jaw/bite opening
- ▶ Jaw position (footplate): 40° and 90° for Ferris–Smith Kerrison rongeur
- ▶ Jaw position: straight, up, down, left, right
- ▶ Footplates: regular and thin footplate
- ▶ Working ends: sharp, sharp semi-serrated (fenestrated), or finely serrated
- ▶ Coated
- ▶ Openable Open Up® versions

To differentiate between the versions, the following additional identifiers are used in conjunction with the catalogue number **REF** RW XXXX:

AU	Punches with pin ejector (AU)
B	TiAlN Ebony (B) coating
O	Open Up® version (openable)
BAJ	Bayonet version
DFP	Punches with thin footplate
E	Punches with external guide (bar ejector)
SK	Inserted part (handle) = “angled”
GEB	Curved version

2.1 Intended use

The rongeurs are used exclusively for removing bone, cartilage, and tissue on the skull and spine.

WARNING

Ferris-Smith Kerrison rongeurs (including bayonet) with thin footplate (DFP)

Article group:

REF

RW 5xxx DFP

may only be used to remove thin or porous bone structures and soft tissue. These Kerrison rongeurs can also be identified by the additional identifier:

The specific suitability of the articles for use on soft tissue, porous bone structures/cartilage and hard bones is listed in the table in Section 1 „Scope“.

The intended user groups are:

1. Surgeons/the specialist responsible for performing the surgery.
2. Support personnel in the operating room:
Refers to those persons who support the surgeon and are expected to perform some of the tasks during the procedure, particularly setting up and adjusting to the surgeon's preferences.
This includes the assistant surgeon, the OR nurses, and the clinical support specialists.

The intended surgical environment includes operating rooms in hospitals, clinics, and outpatient medical centers. The rongeurs are reusable. Their reuse has been validated for 220 applications/reprocessing cycles.

2.2 Indication(s)

The rongeurs are standard instruments for surgical use for neurosurgery and laminectomy.

The treating physician is responsible for selecting the specific rongeur and its surgical use.

2.3 Contraindication(s)

Specific contraindications are not known.

Any complications such as

- local infection
- nerve damage
- failure of the application due to an inadequate healing phase prior to load bearing

are, as a rule, not directly related to the use of the surgical instrument.

3 Safety

3.1 Information about warnings

Warnings refer to the action steps and warn of hazards that may arise when executing the action steps. Warnings precede the action steps.

3.1.1 Warning symbols

The following warning symbols and signal words identify warnings and define the possible consequences of failing to comply with the warnings.

Warning symbol	Signal word	Description	Consequence of noncompliance
	WARNING	Identifies a potential imminent hazard	Mild or moderate injuries may result
	CAUTION	Identifies potential imminent damage to property	The device may be damaged.

3.2 General safety instructions

- ▶ Carefully read the instructions for use
- ▶ Keep instructions for use accessible for users and medical device reprocessing technician
- ▶ Only trained personnel may use the rongeurs.
- ▶ The rongeurs are supplied nonsterile and must undergo the full reprocessing procedure after the transport packaging has been removed.

3.3 Checking before use

WARNING

To avoid a risk of injury and/or malfunction, inspect the rongeur before each use.

- ▶ Check the rongeur for loose, bent, broken, cracked, worn, or broken-off parts.
- ▶ Check that screws are firmly in place.
- ▶ Carry out a functional test.

Do not use instruments with jaw sections with signs of plastic deformation or that do not conform to one another or do not close fully.

3.4 Safety during use

WARNING

Article groups:

REF

All

If excessive forces are applied, unintended injuries to the tissue or bone may cause adverse effects or even break the rongeurs.

Careful use of the instruments is essential.

WARNING

Article groups:

REF

RW 20xx

RW 3xxx

RW 4xxx

Damage or destruction of the jaw sections due to overloading!

- ▶ Only use the rongeurs correctly as described.
- ▶ Do not cut or lever out bone or adjacent bone structures with rongeurs that are only indicated for use on porous bone structures/cartilage and soft tissue.
- ▶ Do not use rongeurs with blunt edges.

- ▶ Avoid overloading the jaw sections by tearing off or holding and levering tissue and avoid applying lateral loads by twisting.
- ▶ Bite volume: The cavity (recess) of the jaw sections determines the bite volume:
The bite volume must therefore be less than the recess on the selected model.
Continuing to use the rongeur despite visible overloading can break the jaw section joint.

WARNING

Article groups:

REF

RW 5xxx

Damage to or destruction of the footplate/cutting area due to overloading!

- ▶ Select the footplate based on the type of application:

Thin footplate:

Instruments with a thin footplate (DFP) may only be used to remove thin or porous bone structures and soft tissue.

These instruments can be identified by the “DFP” suffix added to the article number and the marking “FOR DELICATE USE ONLY”

Regular footplate:

Instruments with a regular footplate are also suitable for removing hard bone and tough tissue.

- ▶ Punch/cutting volume:

The cavity (recess) on the footplate and the sliding shaft limits the punch volume.

The punch/cutting volume must therefore be less than the recess on the selected model.

Bending of the sliding shaft is a sign of overloading:
Stop the punch/cutting process immediately.

Solutions:

Use a rongeur with a wider opening or reduce the punch/cutting volume.

- ▶ Movements while executing punches/cuts:
Avoid overstressing by twisting or levering (particularly with a thin footplate).

If these instructions are not followed, the resulting overloading can cause the footplate and/or the cutting area to bend or break.

CAUTION

Article groups:

REF

All

Damage caused by foreign materials

- ▶ Do not use the rongeurs to hold wires/screws/implants/needles or the like: This may cause dents, ejected material, and material loss in the cutting area.
Do not use rongeurs with damage in the cutting area.

WARNING

Article groups:

REF

Open Up® version (openable)

All with the additional article number identifier „O“

Push button

- ▶ The Open Up® models feature a push button. This is used to place the instruments in the opened-out cleaning position. During surgical use, this push button must not be activated.
- ▶ To prevent inadvertent unlocking, the instrument can only be unlocked when the push button is pushed down fully and the inserted part is pulled (two-hand operation).

If the instrument is nevertheless placed in the opened-out cleaning position during use, the instrument can be locked again: See the Appendix „Handling of the Open Up® System“.

3.5 Reportable incidents

For serious incidents that occur in connection with our medical device, the user is obliged to report the incident to the manufacturer and the competent authorities of the Member State in which the incident occurred.

Please send reports for us to the following email address: reklamation@wenzler.de.

4 Reprocessing procedure

4.1 General instructions

The reprocessing procedure described here has been validated.

The parameters and reprocessing chemicals used in this validation are stated for the reprocessing procedure described here.

4.2 Maximum number of reprocessing cycles

Frequent reprocessing has minimal effects on these instruments.

Their reuse has been validated for 220 applications/reprocessing cycles.

CAUTION

The lifetime of the cutting rongeurs depends primarily on the cutting process being performed carefully by a specialist who correctly uses the rongeurs.

Damage and signs of wear may appear that reduce the service life of the instrument to less than the possible 220 applications/reprocessing cycles.

4.3 Initial treatment at the site of use

If no additional pretreatment is performed, transfer the devices to reprocessing within one hour. Otherwise, dried-on, fixed surgical residues could make cleaning difficult or impossible.

We recommend removing coarse contamination from the instruments immediately after use using a lint-free disposable wipe/paper towel.

CAUTION

- ▶ With wet disposal, do not use fixing agents or warm water (> 40 °C) because they fix residues and may affect the cleaning outcome.
- ▶ The rongeurs are made of stainless steel, and corrosion, pitting, and stress corrosion may occur if the rongeurs are exposed to chlorine or surgical residues containing chlorides (medicinal products) or immersed in saline solution (NaCl).
- ▶ Under no circumstances should oxidizing agents such as H₂O₂ (hydrogen peroxide) be used for black-coated rongeurs (with article number identifier “B”) or gold coatings (e.g., springs). Oxidizing agents can negatively affect the coatings: The titanium component can be leached out by H₂O₂.

4.4 Precleaning

Article groups:

REF

All except those with the additional article number identifier “O”

- ▶ Place the rongeurs in 0.5 % neodisher MediClean forte for 5 minutes and then brush with a soft nylon brush for 2 minutes until the instruments are visibly clean. Move the rongeurs while doing so: Open and close them.
- ▶ Then treat with ultrasound in the cleaning solution for 10 minutes
- ▶ and then rinse under running tap water for 1 minute. Move the rongeurs while doing so: Open and close them.

For the article groups:

REF

Open Up® version (openable)

All with the additional article number identifier “O”

the following procedure is sufficient:

- ▶ Brush the rongeurs with a soft nylon brush for 1 minute while holding under cold, running tap water. While doing so, also place the instrument in the opened-out cleaning position using the push button (see the Appendix “Handling of the Open Up® System”) and treat the inside surfaces.

CAUTION

REF

RW 3xxx O (B)

RW 4xxx O (B)

Avoid applying excessive force to the joint areas of the components while the instrument is in the cleaning position.

- ▶ Avoid applying lateral forces in the marked area:

- ▶ Do not push the sliding shaft past the maximum possible opening.
The area of the pin connection – between the jaw section and the sliding shaft – could be too heavily stressed as a result.

4.5 Automated cleaning and disinfection

Equipment: WD as per the EN ISO 15883 series of standards

Article groups:

REF

Open Up® version (openable)

- ▶ Rongeurs with the Open Up® design (openable) must be cleaned in a washer-disinfector in the opened-out cleaning position

Validated procedure		
Step	Process	
1	Pre-rinse	
	Rinse medium:	Cold mains water
	Exposure time:	1 minute
2	Cleaning	
	Rinse medium:	Cold mains water
	Cleaning temperature:	55 °C
	Cleaning agent:	neodisher®MediClean forte / Manufacturer: Dr. WEIGERT
	Concentration:	At least 0.40 %
	Exposure time:	10 minutes
3	Neutralization (rinsing)	
	Rinse medium:	Cold mains water

Validated procedure		
Step	Process	
	Neutralization agent:	neodisher® N/Manufacturer: Dr. WEIGERT
	Concentration:	0.10 %
	Exposure time:	1 minute
4	Post-rinse	
	Rinse medium:	Warm DI water
	Exposure time:	1 minute
6	Thermal disinfection	
	Disinfection temperature:	At least 90 °C
	Disinfection time:	5 minutes
7	Drying	
	Drying temperature	100 °C
	Drying time	10 minutes

4.6 Inspection

Before the inspection, allow the instruments to cool to room temperature

4.6.1 Visual inspection

After cleaning/disinfection, the instruments must be macroscopically clean, that is, free of visible contamination and residues. The check is performed visually, and critical areas are:

- ▶ sliding surfaces, guides, internal areas, jaw/cutting area, ejector area, push button mechanism, and spring area.
- ▶ If any contamination is detected, the cleaning and disinfection process must be repeated.
- ▶ Check instruments for damage:
 - Cracks, broken-off parts, bent parts
 - Sharp areas, dents, burrs, rough areas on the surface
 - Damage that could damage surgical gloves
 - Check instruments for faded or incomplete labeling

Do not use damaged instruments.

Immediately separate out any instruments that are damaged.

4.6.2 Functional inspection

CAUTION

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REF


Article groups:

All

Damage due to metal wear on the instruments as a result of a lack of lubricant or instrument oil

Before the functional inspection, the movable components must be treated with a lubricant and instrument oil specifically for medical and surgical instruments.

▶ Lubrication points are identified by the  symbol:

For the article groups **REF** Open Up[®] version (openable) with additional article number identifier “O”

Lock an instrument again that is in the opened-out cleaning position; see “Handling of the Open Up[®] System”

- ▶ Check instrument for smooth action and function
- ▶ Check instrument for loose screws

4.7 Packaging

The rongeurs must be placed in a suitable sterile barrier system. The sterile barrier system must meet the following criteria:

- ▶ compliance with the requirements for EN ISO 11607
- ▶ suitable for steam sterilization (steam permeability)

Sterilization accessories and sterilization packaging must be appropriate for the size of the rongeurs as well as the sterilization procedure used.

4.8 Steam sterilization

Steam is the only permissible sterilization medium. The steam must be able to come into contact with all parts of the devices to ensure complete sterilization.

Sterilize the devices in accordance with standard methods for steam sterilization used in hospitals.

Validated sterilization process:

Pre-vacuum	3 times
Sterilization temperature	134 °C
Sterilization time	3 minutes
Drying time	20 minutes

The parameters listed here were selected to cover as broad as possible a range of parameters that are used internationally in sterile processing departments.

Other parameters with a longer holding time and/or higher temperature are therefore also covered.

CAUTION

At the end of the steam sterilization procedure, allow the devices to cool to room temperature (approx. 20 °C) before opening the sterilization packaging.

Accelerated cooling is not permitted and can lead to damage or impaired function.

5 Storage

5.1 Sterile, packaged instruments

The end user is responsible for ensuring that the devices are stored after sterilization so that the sterility of the instruments is maintained until use.

5.2 Nonsterile instruments

The storage period of nonsterile devices is not limited. The devices are manufactured from nondegradable material so that a loss of stability is not expected when the instruments are stored under the recommended conditions.

We recommend storing the devices in a clean and dry environment.

Proper use of the rongeurs is subject to signs of wear and tear. The cutting area in particular is placed under load with every use.

- ▶ Note Section 4.2 “Maximum number of reprocessing cycles” here.

For signs of wear and tear, we offer a technical service: A feasibility check determines whether the extent of the defect/damage permits a repair.

For such returns, use our returns form to make a repair request. The form is available on our website:
www.wenzler.de/service.

You can also use the returns form if you are returning instruments for other reasons.

- ▶ Before returning instruments, used rongeurs must undergo the full reprocessing procedure.
- ▶ Note Section 3.5 “Reportable incidents”.

CAUTION/WARNING

Modifications void the product authorization and warranty claims

Repairs that are not carried out by the manufacturer void the product authorization and warranty claims, therefore:

- ▶ Do not make any modifications
- ▶ Servicing only by Raimund Wenzler GmbH

6.1 Manufacturer/Service address:

Raimund Wenzler GmbH
Instrumenten-Manufaktur
Steigäcker 22
78582 Balgheim / Germany
+49(0)7424 95827-30

7 Disposal

The rongeurs are made of stainless steel and can be recycled.

WARNING

Contamination / Hygiene

- ▶ Before disposal, the rongeurs must undergo the full reprocessing procedure.

WARNING

Danger of injury

Pointed and sharp parts of the rongeurs can cause injuries


- ▶ When disposing of or recycling the instruments, ensure that the packaging prevents injuries.

Symbols used for identification

REF	UDI	GTIN	LOT	SN	MD	QTY
Catalogue number	Unique Device Identifier	Global Trade Item Number	Batch code	Serial number	Medical Device	Quantity
						
Manufacturer	Distributor	Consult instructions for use	Non sterile	Caution	Lubrication points	CE marking
						
Thin footplate	Regular footplate	Sharp jaw section without serration	Jaw section, semi-serrated	Jaw section, finely serrated		

9.1 Trademarks of Raimund Wenzler GmbH

9.2 Trademarks of third-party manufacturers

Trademark	Trademark owner
neodisher® N	Chemische Fabrik Dr. Weigert GmbH & Co. KG
neodisher®MediClean forte	Chemische Fabrik Dr. Weigert GmbH & Co. KG

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- ▶ Laat onderhoud alleen uitvoeren door Raimund Wenzler GmbH
9.1 Merken van Raimund Wenzler GmbH